What Can Clinicians and Patients Expect from Healthpath Gut Health Testing?

0 /

The Healthpath Gut Test shows you what's going on in your gut. By looking at imbalances in bacteria, yeasts, parasites and other intestinal health biomarkers, you find out what's contributing to your symptoms. You also receive targeted diet, supplement and lifestyle recommendations to help you take back control.

The biomarkers provide clinical information on three key areas:



1 | Digestion/Absorption

- pH
- Pancreatic elastase
- Zonulin

2 | Immune activity/Inflammation

- Calprotectin
- Haemoglobin
- Secretory IgA
- H. Pylori
- Archaea/methanogens
- E. Coli, Lactobacillus species, Enterococcus species
- Akkermansia muciniphila, Faecalibacterium prausnitzii



3 | Gut Microbiome/Mycobiome

- Microbiome diversity
- Enterotype
- Dysbiosis index
- Actinobacteria
- Bacteroidetes
- Firmicutes
- Proteobacteria
- Fusobacteria
- Verrucomicrobia
- Hydrogen-sulphide production
- Oxalate-degrading bacteria
- Yeasts/moulds
- Parasites
- Helminths

Clinical Advantages of The Healthpath Gut Health Test qPCR Technology

This new method of analysis allows for a single sample. This makes the process easier for everyone, and it's particularly helpful for children and those struggling with diarrhoea or constipation.



Stool properties		
Colour	\checkmark	\checkmark
Consistency	\checkmark	\checkmark
рН	\checkmark	\checkmark

Diversity:

Your diversity is key, which is why our microbiome analysis covers hundreds of parameters. High bacterial diversity is known to protect against intestinal infections. But low bacterial diversity is common, especially in disease states or after a course of antibiotics. When diversity is low, opportunistic bacteria like pathogens, fungi and viruses can proliferate.

Rather than focusing on individual species, it's more important to investigate how the different bacteria interact. Together, they're responsible for a host of intestinal functions.

Biodiversity		
Diversity	\checkmark	\checkmark
Dysbiosis index	\checkmark	\checkmark

There are four large phyla (groups) of bacteria: Bacteroidetes, Firmicutes, Actinobacteria and Proteobacteria. We also report on two smaller, clinically relevant phyla: Verrucomicrobia and Fusobacteria.

Bacterial distribution

Actinobacteria	\checkmark	\checkmark
Bacteroidetes	\checkmark	\checkmark
Firmicutes	\checkmark	\checkmark
Fusobacteria	\checkmark	\checkmark
Proteobacteria	\checkmark	\checkmark
Verrucomicrobia	\checkmark	\checkmark
Other	\checkmark	\checkmark
Firmicutes/ Bacteroidetes Ratio	\checkmark	\checkmark

Enterotype:

Recent research suggests there are three different types of gut microbiomes, known as 'enterotypes'. Not only do the different enterotypes influence the absorption of minerals, but they also have different metabolic properties.

Advanced

Gut

Health

Test Pro

Gut

Health

Test

Enterotype 1 has high levels of Bacteroides species, which use fat and protein effectively. Enterotype 2 has a strong Prevotella population, which is better at metabolising carbohydrates. Enterotype 3 is the rarest enterotype. It has high levels of Ruminococcus flora, though we don't yet know which macronutrients it prefers.

Enterotypes aren't affected by a person's age or gender and they remain stable for years. They can be influenced, however, by a long-term change of diet and by taking prebiotics.

Enterotype

1, 2 or 3	\checkmark	\checkmark

Actinobacteria

Bifidobacteria		\checkmark	\checkmark		
Equol-producing bacteria		\checkmark	\checkmark		
	Adlercreutzia species		\checkmark		
	Eggerthella lenta		\checkmark		
	Slackia species		\checkmark		

Bacteroidetes			
Ba	cteroides	\checkmark	\checkmark
Prevotella		\checkmark	\checkmark
	Prevotella copri	\checkmark	\checkmark



Firmicutes:

Butyrate is a short-chain fatty acid that's produced by bacteria in the colon. It's quickly absorbed by the intestinal mucosa, which means the only reliable way to measure it is to look at the number of butyrate-producing bacteria.

Firmicutes bacteria are key butyrate producers. One of these, Faecalibacterium prausnitzii, typically makes up 5–15% of human intestinal bacteria. This important butyrate-producing species has anti-inflammatory properties—so much so that an absence of Faecalibacterium prausnitzii typically correlates with higher levels of inflammation.

Firmicutes

Butyrate-producing bacteria		\checkmark	\checkmark
	Faecalibacterium prausnitzii	\checkmark	\checkmark
	Eubacterium rectale	\checkmark	\checkmark
	Eubacterium hallii	\checkmark	\checkmark
	Roseburia species	\checkmark	\checkmark
	Ruminococcus species	\checkmark	\checkmark
	Coprococcus	\checkmark	\checkmark
	Butyrivibrio species		\checkmark
	Cl. butyricum		\checkmark
	Total bacterial count	\checkmark	\checkmark
Clostridia		\checkmark	\checkmark
	Clostridia total bacterial count	\checkmark	\checkmark
	Clostridia cluster 1	\checkmark	\checkmark
	Clostridia histolytium		\checkmark
	Clostridium perfringens		\checkmark
	Clostridium sporenges		\checkmark
Other			\checkmark
	Christensenellaceae		\checkmark
	Dialister invisus		\checkmark

Fusobacteria		
Fusobacterium species	\checkmark	\checkmark

Verrucomicrobia		
Akkermansia muciniphila	\checkmark	\checkmark

Gut	Advance
Health	Gut
Test	Health
Test	Test Pr

Proteobacteria

Potentially pathogenic bacteria		ntially pathogenic bacteria	\checkmark	\checkmark
	Haemophilus		\checkmark	\checkmark
	Acinetobacter		\checkmark	\checkmark
	Escherichia coli biovare		\checkmark	\checkmark
	Pro	oteus species	\checkmark	\checkmark
		Proteus mirabilis		\checkmark
	Kle	ebsiella species	\checkmark	\checkmark
		Klebsiella pneumoniae		\checkmark
	Enterobacter species		\checkmark	\checkmark
	Se	rratia species	\checkmark	\checkmark
	На	fnia species	\checkmark	\checkmark
	Mc	organella species	\checkmark	\checkmark
	Campylobacter species			\checkmark
	Pro	ovidencia species		\checkmark
	Cit	robacter species		\checkmark
Hi	Histamine-producing bacteria		\checkmark	\checkmark
H2S production		production	\checkmark	\checkmark

Hydrogen-sulphide production:

Bacterial metabolism isn't always a good thing. Some bacteria reduce sulphate to create hydrogen sulphide—a toxic metabolic by-product that can damage the gut lining. The species Bilophila wadsworthii, Desulfomonas pigra and Desulfovibrio piger are thought to be potent hydrogen-sulphide developers.

	Sulphate-reducing bacteria	\checkmark	\checkmark
	Desulfovibrio piger		\checkmark
	Desulfomonas pigra		\checkmark
	Bilophila wadsworthii		\checkmark
Oxalate-degrading bacteria			\checkmark
	Oxalobacter formigenes		\checkmark

Archaea:

Archaea have been overlooked in microbiome studies until recently. New research suggests that 1) archaea are part of the microbiome in plants, animals and humans, 2) they form biofilms and 3) they interact with the human immune system. Some archaea are also methanogens, which may play a role in chronic constipation.

 \checkmark

 \checkmark

Archaea

Methanobrevibacter



Immunogenically effective bacteria		
Escherichia coli	\checkmark	\checkmark
Enterococcus species	\checkmark	\checkmark
Lactobacillus species	\checkmark	\checkmark

Mucin production/mucosal barrier:

A healthy colon has a protective mucous layer. If this layer is damaged—or only small amounts of mucous are produced—pathogens, pollutants and allergens can come into direct contact with the mucosa. This leads to inflammation.

The bacterium Akkermansia muciniphila is important because it encourages goblet cells to produce this protective mucous. Parts of this mucous also provide a special type of carbohydrate called oligosaccharides, which feed the bacteria that make gut-healing butyrate. With the right bacteria, it becomes a virtuous circle!

Mucin production/ mucosal barrier

Akkermansia muciniphila	\checkmark	\checkmark
Faecalibacterium prausnitzii	\checkmark	\checkmark

Yeasts/moulds		
Candida albicans	\checkmark	\checkmark
Candida species	\checkmark	\checkmark
Geotrichum candidum	\checkmark	\checkmark
Moulds	\checkmark	\checkmark

Functional markers

Calprotectin	\checkmark	\checkmark
Haemoglobin in faeces immunologically	\checkmark	\checkmark
Secretory IgA	\checkmark	\checkmark
Pancreatic elastase	\checkmark	\checkmark
Zonulin		\checkmark

Gut Health Test Advanced Gut Health Test Pro

Parasites:

The Multiplex Real-time PCR (Multiplex quantitative real-time PCR) is a faster and more effective method for detecting parasites. This new test:

- provides reliable analysis, even with minimal attack
- gives no false positives with non-pathogens
- can be sent out with regular mail
- gives reliable results in symptom-free patients and also after treatment

Parasites

Pathobionts		\checkmark	\checkmark
	Blastocystis hominis	\checkmark	\checkmark
	Dientamoeba fragilis	\checkmark	\checkmark
	Helicobacter AG	\checkmark	\checkmark
Pa	athogenic intestinal protozoa	\checkmark	\checkmark
	Giardia lamblia	\checkmark	\checkmark
	Entamoeba histolytica	\checkmark	\checkmark
	Cryptosporidium species	\checkmark	\checkmark
	Cyclospora cayetanensis	\checkmark	\checkmark
Helminths COMING SOON			\checkmark
	Taenia species		\checkmark
	Taenia solium		\checkmark
	Taenia saginata		\checkmark
	Ascaris species		\checkmark
	Enterobius vermicularis		\checkmark
	Ancylostoma species		\checkmark
	Ancylostoma duodenale		\checkmark
	Hymenolepsis species		\checkmark
	Hymenolepsis nana		\checkmark
	Hymenolepsis diminuta		\checkmark
	Trichuris trichiura		\checkmark
	Necator americanus		\checkmark
	Strongyloides species		\checkmark
	Strongyloides stercoralis		\checkmark
	Microsporidia		\checkmark
	Enterocytozoon species		\checkmark
	Encephalitozoon species		\checkmark